



## Product description

Order number	93.854.123
Product description	Block Rack D12, for 50 collection tubes, Ø opening: 12 mm, format: 5 x 10, green, material: PP

## Product characteristics

Number of wells	50
Storage capacity	50
Format	5 x 10
Autoclavable (max. °C)	121 °C

## Size

Diameter of opening	12 mm
Width of product	80 mm
Height of product	45 mm
Length of product	152 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	green

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	1
Type of smallest subpackaging	-
Piece(s) / inner box	1
Piece (s) / outer case	20
Piece(s) / pallet	1560
Depth of case	428 mm
Width of case	323 mm
Height of case	130 mm
Case volume	0.018 cbm
Weight of product	0.1214 kg
Weight of case	2.8 kg
EAN of inner box	4038917122665
EAN case	4038917018500

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

