



Product Specification
CryoRack 40, PP, format: 10 x 4, suitable for
CryoPure tubes

Page 1



Product description

Order number	93.856.040
Product description	CryoRack 40, material: PP, format: 10 x 4, (LxWxH): 100 x 210 x 22 mm, for 40 collection tubes, suitable for CryoPure tubes, 1 piece(s)/bag

Product characteristics

Graduation	no
Suitable for	CryoPure tubes
Storage capacity	40
Format	10 x 4

Size

Width of product	210 mm
Height of product	22 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest
subpackaging bag
Piece(s) / inner box 1
Piece (s) / outer case 10
Piece(s) / pallet 1650
Depth of case 228 mm
Width of case 218 mm
Height of case 155 mm
Case volume 0.0077 cbm
Weight of product 0.0679 kg
Weight of case 0.9 kg
EAN of inner box 4038917122788
EAN case 4038917077293

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

