



### **Product description**

Order number	93.893.100
Product description	Rack, material: PC, transparent, format: 10 x 2, (LxWxH): 327 x 72 x 60 mm, for 20 collection tubes, suitable for tubes 100 x 21.5 mm, 1 piece(s)/case

### **Product characteristics**

Suitable for	tubes 100 x 21.5 mm
Storage capacity	20
Format	10 x 2
Autoclavable (max. °C)	121 °C

### **Size**

Width of product	72 mm
Height of product	60 mm
Length of product	327 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polycarbonate (PC)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	40
Piece(s) / pallet	800
Depth of case	584 mm
Width of case	387 mm
Height of case	360 mm
Case volume	0.0814 cbm
Weight of product	0.2524 kg
Weight of case	11.57 kg
EAN of inner box	4038917123136
EAN case	4038917362375

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

