




Product Specification
Rack, styrofoam, format: 20 x 5, suitable for
S-Monovette® Ø 11 mm

Page 1



Product description

| | |
|---------------------|--|
| Order number | 95.1439 |
| Product description | Rack, material: styrofoam, white, format: 20 x 5, (LxWxH): 338 x 90 x 50 mm, suitable for S-Monovette® Ø 11 mm, 1 piece(s)/blister |

Product characteristics

| | |
|--------------|----------------------|
| Suitable for | S-Monovette® Ø 11 mm |
| Format | 20 x 5 |

Size

| | |
|-------------------|--------|
| Width of product | 90 mm |
| Height of product | 50 mm |
| Length of product | 338 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Rack, styrofoam, format: 20 x 5, suitable for
S-Monovette® Ø 11 mm

Page 2

Material & colours

| | |
|-------------------|-----------|
| Product material | styrofoam |
| Colour of product | white |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | blister |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 40 |
| Piece(s) / pallet | 640 |
| Depth of case | 498 mm |
| Width of case | 428 mm |
| Height of case | 380 mm |
| Case volume | 0.081 cbm |
| Weight of product | 0.0235 kg |
| Weight of case | 1.7 kg |
| EAN of inner box | 4038917086646 |
| EAN case | 4038917005173 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT