



Product description

| | |
|---------------------|---|
| Order number | 95.64.251 |
| Product description | Rack, material: styrofoam, white, format: 10 x 5, (LxWxH): 205 x 127 x 46 mm, for 50 collection tubes, suitable for 15 ml tubes, 1 piece(s)/blister |

Product characteristics

| | |
|------------------|-------------|
| Suitable for | 15 ml tubes |
| Storage capacity | 50 |
| Format | 10 x 5 |

Size

| | |
|-------------------|--------|
| Width of product | 127 mm |
| Height of product | 46 mm |
| Length of product | 205 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|-----------|
| Product material | styrofoam |
| Colour of product | white |
| Closure material | styrofoam |
| Colour of cap | white |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 21 |
| Type of smallest subpackaging | blister |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 21 |
| Piece(s) / pallet | 882 |
| Depth of case | 388 mm |
| Width of case | 378 mm |
| Height of case | 235 mm |
| Case volume | 0.0345 cbm |
| Weight of product | 0.0271 kg |
| Weight of case | 0.58 kg |
| EAN of inner box | 4038917095594 |
| EAN case | 4038917007771 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

