



### **Product description**

Order number	95.64.253
Product description	Rack, lid enclosed, material: styrofoam, white, format: 5 x 4, (LxWxH): 100 x 85 x 65 mm, for 20 collection tubes, suitable for tubes Ø 10.8 mm, 1 piece(s)/blister

### **Product characteristics**

Suitable for	tubes Ø 10.8 mm
Storage capacity	20
Format	5 x 4

### **Size**

Width of product	85 mm
Height of product	65 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	styrofoam
Colour of product	white
Closure material	styrofoam
Colour of cap	white

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	112
Piece(s) / pallet	1792
Depth of case	498 mm
Width of case	428 mm
Height of case	380 mm
Case volume	0.081 cbm
Weight of product	0.0125 kg
Weight of case	1.96 kg
EAN of inner box	4038917095600
EAN case	4038917007788

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

