




Product Specification
Cryobox, 132 x 132 x 53 mm, format: 10 x 10, for
100 collection tubes

Page 1



Product description

Order number	93.874.410
Product description	Cryobox, with numerical coding at each aperture, for low-temperature storage, material: PC, yellow, slip-on lid with ventilation function, cap: transparent, (LxWxH): 132 x 132 x 53 mm, format: 10 x 10, for 100 collection tubes, for CryoPure tubes 1.2?2.0 ml internal thread, 5 piece(s)/bag

Product characteristics

Graduation	no
Cap	slip-on lid with ventilation function
Closure type	slip-on lid
Design	with numerical coding at each aperture
Suitable for	CryoPure tubes 1.2?2.0 ml internal thread
Storage capacity	100
Format	10 x 10

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

Width of product	132 mm
Height of product	53 mm
Length of product	132 mm

Material & colours

Product material	Polycarbonate (PC)
Colour of product	yellow
Colour of cap	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	5
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	20
Piece(s) / pallet	720
Depth of case	135 mm
Width of case	270 mm
Height of case	540 mm
Case volume	0.0197 cbm
Weight of product	0.1593 kg
Weight of case	3.71 kg
EAN of inner box	4038917122917
EAN case	4038917415163

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

