



Product Specification
Storage box, hinged lid, PP, format: 10 x 10, for 100 collection tubes

Page 1



Product description

Order number	93.877.999
Product description	Storage box, hinged lid, material: PP, colour mix, format: 10 x 10, for 100 collection tubes, suitable for vessels with maximum dimensions of 45 x 12 mm

Product characteristics

Closure type	hinged lid
Suitable for	vessels with maximum dimensions of 45 x 12 mm
Storage capacity	100
Format	10 x 10
Autoclavable (max. °C)	121 °C

Size

Width of product	140 mm
Height of product	55 mm
Length of product	150 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	colour mix
Closure material	Polypropylene (PP)
Colour of cap	colour mix

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	5
Piece(s) / inner box	5
Piece (s) / outer case	20
Piece(s) / pallet	720
Case volume	0.026 cbm
Weight of product	0.204 kg
Weight of case	4.09 kg
EAN of inner box	4038917240376
EAN case	4038917240369

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

