



Product description

Order number	93.881
Product description	Storage box, slip-on lid, material: PP, blue, format: 4 x 4, for 16 collection tubes, suitable for 50 ml tubes and tubes up to Ø 30 mm, 2 piece(s)/bag

Product characteristics

Closure type	slip-on lid
Suitable for	50 ml tubes and tubes up to Ø 30 mm
Storage capacity	16
Format	4 x 4
Autoclavable (max. °C)	121 °C

Size

Width of product	144 mm
Height of product	127 mm
Length of product	144 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	blue
Closure material	Polypropylene (PP)
Colour of cap	transparent

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	2
Type of smallest subpackaging	bag
Piece(s) / inner box	2
Piece (s) / outer case	12
Piece(s) / pallet	420
Depth of case	251 mm
Width of case	300 mm
Height of case	452 mm
Case volume	0.034 cbm
Weight of product	0.2228 kg
Weight of case	3.2 kg
EAN of inner box	4038917355223
EAN case	4038917355216

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

