



Product description

| | |
|---------------------|--|
| Order number | 95.064.259 |
| Product description | Box, Styrofoam, 20-place, for micro tubes, 8 pcs./pack |

Product characteristics

| | |
|------------------|-------|
| Storage capacity | 20 |
| Format | 5 x 4 |

Material & colours

| | |
|-------------------|-----------|
| Product material | Styrofoam |
| Colour of product | white |

Purity & certification

| | |
|------------------|----------------------|
| Product category | No marking necessary |
| Batched | Lot number |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

| | |
|------------------------|---------------|
| Minimum order qty. | 16 |
| Piece(s) / inner box | 8 |
| Piece (s) / outer case | 56 |
| Piece(s) / pallet | 1792 |
| Depth of case | 401 mm |
| Width of case | 451 mm |
| Height of case | 280 mm |
| Case volume | 0.0504 cbm |
| Weight of case | 1.36 kg |
| EAN of inner box | 4038917386883 |
| EAN case | 4038917386876 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

