



Product description

Order number	95.64.259
Product description	Rack, material: styrofoam, white, format: 5 x 4, for 20 collection tubes, 8 piece(s)/case

Product characteristics

Storage capacity	20
Format	5 x 4

Material & colours

Product material	styrofoam
Colour of product	white

Purity & certification

Product category	no medical device IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Rack, styrofoam, format: 5 x 4

Page 2

Packaging

Minimum order qty.	56
Type of smallest subpackaging	case
Piece(s) / inner box	8
Piece (s) / outer case	56
Piece(s) / pallet	1792
Case volume	0.0559 cbm
Weight of product	0.0238 kg
Weight of case	1.36 kg
EAN of inner box	4038917386883
EAN case	4038917386876

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

