



Product Specification
**Storage box, slip-on lid, cardboard, format: 4 x 4,
for 16 collection tubes**



Product description

Order number	95.64.916
Product description	Storage box, slip-on lid, material: cardboard, white, format: 4 x 4, for 16 collection tubes, suitable for Tubes up to max. 34 mm Ø, tube height between 89-115 mm, 1 piece(s)/case

Product characteristics

Closure type	slip-on lid
Suitable for	Tubes up to max. 34 mm Ø, tube height between 89-115 mm
Storage capacity	16
Format	4 x 4

Size

Width of product	148 mm
Height of product	116 mm
Length of product	148 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	cardboard
Colour of product	white
Closure material	cardboard
Colour of cap	white

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	30
Piece(s) / pallet	360
Depth of case	495 mm
Width of case	419 mm
Height of case	546 mm
Case volume	0.1134 cbm
Weight of product	0.1337 kg
Weight of case	5.15 kg
EAN of inner box	4038917095679
EAN case	4038917007832

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

