



**Product Specification**  
**Storage box, slip-on lid, cardboard, format: 7 x 7,  
for 49 collection tubes**

**Page 1**



**Product description**

Order number	95.64.950
Product description	Storage box, slip-on lid, material: cardboard, white, format: 7 x 7, for 49 collection tubes, suitable for tubes Ø 16-17 mm

**Product characteristics**

Closure type	slip-on lid
Suitable for	tubes Ø 16-17 mm
Storage capacity	49
Format	7 x 7

**Size**

Width of product	135 mm
Height of product	87 mm
Length of product	135 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Material & colours**

Product material	cardboard
Colour of product	white
Closure material	cardboard
Colour of cap	white

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	60
Piece(s) / inner box	1
Piece (s) / outer case	60
Piece(s) / pallet	960
Case volume	0.1134 cbm
Weight of product	0.1373 kg
Weight of case	8.27 kg
EAN of inner box	4038917319607
EAN case	4038917319591

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

