



Product Specification
Push cap, white, suitable for 2 ml sample tube
73.663



Product description

Order number 65.108
Product description Push cap, white, suitable for 2 ml sample tube 73.663,
1,000 piece(s)/bag

Product characteristics

Cap flat cap
Closure type push cap
Suitable for 2 ml sample tube 73.663

Material & colours

Product material Low Density Polyethylene (LD-PE)
Colour of product white

Purity & certification

Product category no medical device | IVD
Batched no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	10000
Piece(s) / pallet	640000
Depth of case	398 mm
Width of case	298 mm
Height of case	215 mm
Case volume	0.0255 cbm
Weight of product	0.0004 kg
Weight of case	4.1 kg
EAN of inner box	4038917080866
EAN case	4038917000048

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

