



Product Specification
Ventilation cap, natural, suitable for tubes Ø 12 mm

Page 1



Product description

Order number	65.722
Product description	Ventilation cap, natural, suitable for tubes Ø 12 mm, 100 piece(s)/bag

Product characteristics

Cap	Ventilation cap
Closure type	push cap
Suitable for	tubes Ø 12 mm

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

Purity & certification

Product category	no medical device IVD
Latex-free	yes
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Ventilation cap, natural, suitable for tubes Ø 12 mm

Page 2

Packaging

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	5000
Piece(s) / pallet	240000
Depth of case	388 mm
Width of case	378 mm
Height of case	215 mm
Case volume	0.0315 cbm
Weight of product	0.0006 kg
Weight of case	3.71 kg
EAN of inner box	4038917083706
EAN case	4038917002608

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-06

This document was prepared by EDP support and is valid without signature.



SARSTEDT