



**Product Specification**  
**Ventilation cap, natural, suitable for tubes Ø 16 mm**

**Page 1**



**Product description**

Order number	65.725
Product description	Ventilation cap, natural, suitable for tubes Ø 16 mm, 100 piece(s)/bag

**Product characteristics**

Cap	Ventilation cap
Closure type	push cap
Suitable for	tubes Ø 16 mm

**Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

**Purity & certification**

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-05

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Ventilation cap, natural, suitable for tubes Ø 16 mm**

**Page 2**

**Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	5000
Piece(s) / pallet	120000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.001 kg
Weight of case	5.88 kg
EAN of inner box	4038917083737
EAN case	4038917002639

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-05

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**