



**Product Specification**  
**Tube, 5 ml, (LxØ): 75 x 12 mm, PS**

**Page 1**



**Product description**

Order number	55.1579.002
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PS, round base, transparent, push cap, suitable for flow cytometer like FACSCanto, FACS Aria, FACSCalibur, LRS2, FC500, sterile, 50 piece(s)/bag

**Product characteristics**

Graduation	no
Closure type	push cap
Suitable for	flow cytometer like FACSCanto, FACS Aria, FACSCalibur, LRS2, FC500
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	4000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tube, 5 ml, (LxØ): 75 x 12 mm, PS**

**Page 2**

**Size**

Diameter	12 mm
Length of product	75 mm

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	1000
Piece(s) / pallet	56000
Depth of case	598 mm
Width of case	398 mm
Height of case	125 mm
Case volume	0.0298 cbm
Weight of product	0.0026 kg
Weight of case	3.18 kg
EAN of inner box	4038917162746
EAN case	4038917162739

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**