



Product description

Order number	55.461
Product description	Tube, working volume: 14 ml, (LxØ): 119 x 15.5 mm, material: PS, round base, transparent, push cap, 500 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	14 ml
Centrifugation max (RCF)	3500 x g

Size

Diameter	15.5 mm
Length of product	119 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 14 ml, (LxØ): 119 x 15.5 mm, PS

Page 2

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	30000
Depth of case	292 mm
Width of case	464 mm
Height of case	413 mm
Case volume	0.0559 cbm
Weight of product	0.0058 kg
Weight of case	6.29 kg
EAN of inner box	4038917081696
EAN case	4038917000802

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-13

This document was prepared by EDP support and is valid without signature.



SARSTEDT