



### Product description

Order number	55.467.005
Product description	Tube, working volume: 8 ml, (LxØ): 100 x 13 mm, material: PS, round base, transparent, push cap, 500 piece(s)/bag

### Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	8 ml
Centrifugation max (RCF)	4000 x g

### Size

Diameter	13 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tube, 8 ml, (LxØ): 100 x 13 mm, PS**

**Page 2**

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	48000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0041 kg
Weight of case	9.05 kg
EAN of inner box	4038917106825
EAN case	4038917011570

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-21

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**