



Product description

| | |
|---------------------|--|
| Order number | 55.475.030 |
| Product description | Tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PS, round base, transparent, 500 piece(s)/StackPack |

Product characteristics

| | |
|--------------------------|------------|
| Graduation | no |
| Base shape | round base |
| Volume of work | 5 ml |
| Centrifugation max (RCF) | 4000 x g |

Size

| | |
|-------------------|-------|
| Diameter | 13 mm |
| Length of product | 75 mm |

Material & colours

| | |
|-------------------|------------------|
| Product material | Polystyrene (PS) |
| Colour of product | transparent |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 5 ml, (LxØ): 75 x 13 mm, PS

Page 2

Purity & certification

| | |
|---------------------------|--------------------------------|
| Satisfies the requirement | sealed at 0.6 bar/10 min |
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 2000 |
| Type of smallest subpackaging | StackPack |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 112000 |
| Depth of box | 416 mm |
| Width of box | 77 mm |
| Height of box | 193 mm |
| Depth of case | 428 mm |
| Width of case | 383 mm |
| Height of case | 175 mm |
| Case volume | 0.0287 cbm |
| Weight of product | 0.0031 kg |
| Weight of case | 6.89 kg |
| EAN of inner box | 4038917357883 |
| EAN case | 4038917357876 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT