




Product Specification
Tube, 8.5 ml, (LxØ): 75 x 15.7 mm, PS

Page 1



Product description

Order number	55.495
Product description	Tube, working volume: 8.5 ml, (LxØ): 75 x 15.7 mm, material: PS, round base, transparent, push cap, without cap, 1,000 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	8,5 ml
Centrifugation max (RCF)	4000 x g

Size

Diameter	15.7 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 8.5 ml, (LxØ): 75 x 15.7 mm, PS

Page 2

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	60000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.0042 kg
Weight of case	9.01 kg
EAN of inner box	4038917082013
EAN case	4038917001120

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT