



### Product description

|                     |  |
|---------------------|--|
| Order number        | 55.525.061   |
| Product description | Tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PP, round base, transparent, push cap, with graduations, 100 piece(s)/bag |

### Product characteristics

|                          |            |
|--------------------------|------------|
| Graduation               | yes        |
| Closure type             | push cap   |
| Base shape               | round base |
| Volume of work           | 5 ml       |
| Centrifugation max (RCF) | 5000 x g   |
| Autoclavable (max. °C)   | 121 °C     |

### Size

|                   |       |
|-------------------|-------|
| Diameter          | 13 mm |
| Height of product | 75 mm |
| Length of product | 75 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tube, 5 ml, (LxØ): 75 x 13 mm, PP**

**Page 2**

**Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |

**Purity & certification**

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Batched          | yes                           |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 4000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 100           |
| Piece (s) / outer case        | 4000          |
| Piece(s) / pallet             | 48000         |
| Case volume                   | 0.1237 cbm    |
| Weight of product             | 0.001 kg      |
| Weight of case                | 5.85 kg       |
| EAN of inner box              | 4038917276825 |
| EAN case                      | 4038917276818 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**