



**Product Specification**  
**Tube, 5 ml, (LxØ): 75 x 12 mm, PP**

**Page 1**



**Product description**

|                     |   |
|---------------------|---|
| Order number        | 55.526.521  |
| Product description | Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, blue, push cap, 500 piece(s)/StackPack |

**Product characteristics**

|                          |            |
|--------------------------|------------|
| Graduation               | no         |
| Closure type             | push cap   |
| Base shape               | round base |
| Volume of work           | 5 ml       |
| Centrifugation max (RCF) | 5000 x g   |
| Autoclavable (max. °C)   | 121 °C     |

**Size**

|                   |       |
|-------------------|-------|
| Diameter          | 12 mm |
| Height of product | 75 mm |
| Length of product | 75 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tube, 5 ml, (LxØ): 75 x 12 mm, PP**

**Page 2**

**Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | blue               |

**Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | yes                     |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 500           |
| Type of smallest subpackaging | StackPack     |
| Piece(s) / inner box          | 500           |
| Piece (s) / outer case        | 4000          |
| Piece(s) / pallet             | 120000        |
| Case volume                   | 0.046 cbm     |
| Weight of product             | 0.0021 kg     |
| Weight of case                | 10.25 kg      |
| EAN of inner box              | 4038917109758 |
| EAN case                      | 4038917031950 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**