



### **Product description**

Order number	55.526.524
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, orange, push cap, 500 piece(s)/StackPack

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g
Autoclavable (max. °C)	121 °C

### **Size**

Diameter	12 mm
Height of product	75 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tube, 5 ml, (LxØ): 75 x 12 mm, PP**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	orange

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	StackPack
Piece(s) / inner box	500
Piece (s) / outer case	4000
Piece(s) / pallet	120000
Case volume	0.046 cbm
Weight of product	0.0021 kg
Weight of case	10.25 kg
EAN of inner box	4038917109789
EAN case	4038917032209

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**