



Product Specification
Tube, 4 ml, (LxØ): 50 x 14 mm, PS

Page 1



Product description

Order number	55.655
Product description	Tube, working volume: 4 ml, (LxØ): 50 x 14 mm, material: PS, round base, transparent, push cap, 1,000 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	4 ml
Centrifugation max (RCF)	5000 x g

Size

Diameter	14 mm
Length of product	50 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 4 ml, (LxØ): 50 x 14 mm, PS

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	3000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	3000
Piece(s) / pallet	108000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.002 kg
Weight of case	6.78 kg
EAN of inner box	4038917083157
EAN case	4038917002127

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

