



Product description

Order number	57.512.001
Product description	Tube, working volume: 4.5 ml, (LxØ): 75 x 12 mm, material: PP, conical base, transparent, push cap, without cap, 500 piece(s)/StackPack

Product characteristics

Graduation	no
Closure type	push cap
Base shape	conical base
Volume of work	4,5 ml
Centrifugation max (RCF)	5000 x g
Autoclavable (max. °C)	121 °C

Size

Diameter	12 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	50000
Type of smallest subpackaging	StackPack
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	120000
Depth of box	387 mm
Width of box	77 mm
Height of box	176 mm
Depth of case	398 mm
Width of case	358 mm
Height of case	175 mm
Case volume	0.0249 cbm
Weight of product	0.0021 kg
Weight of case	4.88 kg
EAN of inner box	4038917108799
EAN case	4038917072489

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

