



## Product description

Order number	57.527.300
Product description	Tube, working volume: 12 ml, (LxØ): 110 x 17 mm, material: PP, conical base, transparent, push cap, with graduations, 500 piece(s)/bag

## Product characteristics

Graduation	yes
Closure type	push cap
Base shape	conical base
Volume of work	12 ml
Centrifugation max (RCF)	3000 x g
Autoclavable (max. °C)	121 °C

## Size

Diameter	17 mm
Height of product	110 mm
Length of product	110 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tube, 12 ml, (LxØ): 110 x 17 mm, PP**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Case volume	0.0492 cbm
Weight of product	0.005 kg
Weight of case	5.48 kg
EAN of inner box	4038917367431
EAN case	4038917367424

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-10

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**