



Product Specification
Tube, 7 ml, (LxØ): 50 x 16 mm, PS

Page 1



Product description

Order number	58.485.300
Product description	Tube, working volume: 7 ml, (LxØ): 50 x 16 mm, material: PS, flat base, transparent, push cap, 1,000 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	flat base
Volume of work	7 ml

Size

Diameter	16 mm
Height of product	50 mm
Length of product	50 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Tube, 7 ml, (LxØ): 50 x 16 mm, PS

Page 2

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	72000
Case volume	0.0492 cbm
Weight of product	0.0028 kg
Weight of case	6.07 kg
EAN of inner box	4038917367462
EAN case	4038917367455

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT