



Product description

| | |
|---------------------|---|
| Order number | 58.485 |
| Product description | Tube, working volume: 7 ml, (LxØ): 50 x 16 mm, material: PS, flat base, transparent, push cap, 1,000 piece(s)/bag |

Product characteristics

| | |
|--------------------------|-----------|
| Graduation | no |
| Closure type | push cap |
| Base shape | flat base |
| Volume of work | 7 ml |
| Centrifugation max (RCF) | 2500 x g |

Size

| | |
|-------------------|-------|
| Diameter | 16 mm |
| Length of product | 50 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 7 ml, (LxØ): 50 x 16 mm, PS

Material & colours

| | |
|-------------------|------------------|
| Product material | Polystyrene (PS) |
| Colour of product | transparent |

Purity & certification

| | |
|---------------------------|--------------------------|
| Satisfies the requirement | sealed at 0.6 bar/10 min |
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 2000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1000 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 72000 |
| Depth of case | 398 mm |
| Width of case | 398 mm |
| Height of case | 300 mm |
| Case volume | 0.0475 cbm |
| Weight of product | 0.0028 kg |
| Weight of case | 6.34 kg |
| EAN of inner box | 4038917081924 |
| EAN case | 4038917001038 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

