



Product description

Order number	58.505.300
Product description	Tube, working volume: 12 ml, (LxØ): 95 x 16.5 mm, material: PP, flat base, transparent, push cap, 500 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	flat base
Volume of work	12 ml
Autoclavable (max. °C)	121 °C

Size

Diameter	16.5 mm
Height of product	95 mm
Length of product	95 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 12 ml, (LxØ): 95 x 16.5 mm, PP

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Case volume	0.0492 cbm
Weight of product	0.0049 kg
Weight of case	5.4 kg
EAN of inner box	4038917367523
EAN case	4038917367516

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT