



Product Specification

Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, without cap, 1,000 piece(s)/bag

Page 1



Product description

Order number	60.504.010
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 75 x 13 mm, round base, transparent, material: PP, without cap, 1,000 piece(s)/bag, 2,000 piece(s)/carton

Product characteristics

Graduation	no
Closure type	screw cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g

Size

Diameter	13 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, without cap, 1,000 piece(s)/bag

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Tube material	Polypropylene (PP)

Purity & certification

Satisfies the requirement	IATA
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	72000
Depth of case	398 mm
Width of case	398 mm
Height of case	265 mm
Case volume	0.042 cbm
Weight of product	0.0024 kg
Weight of case	5.43 kg
EAN of inner box	4038917309059
EAN case	4038917531764

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT