



Product description

| | |
|---------------------|---|
| Order number | 60.540.352 |
| Product description | Screw cap tube, working volume: 13 ml, (LxØ): 101 x 16.5 mm, material: PP, round base, transparent, screw cap, blue, cap assembled, sterile, 500 piece(s)/bag |

Product characteristics

| | |
|--------------------------|------------|
| Graduation | no |
| Closure type | screw cap |
| Base shape | round base |
| Volume of work | 13 ml |
| Centrifugation max (RCF) | 4000 x g |

Size

| | |
|-------------------|---------|
| Diameter | 16.5 mm |
| Length of product | 101 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | blue |

Purity & certification

| | |
|---------------------------|-------------------------|
| Satisfies the requirement | IATA |
| Product category | no medical device IVD |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 28000 |
| Case volume | 0.029 cbm |
| Weight of product | 0.006 kg |
| Weight of case | 3.45 kg |
| EAN of inner box | 4038917367646 |
| EAN case | 4038917367639 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

