



Product description

Order number	60.550
Product description	Screw cap tube, working volume: 7 ml, (LxØ): 82 x 13 mm, material: PP, flat base, transparent, screw cap, natural, cap enclosed, 1,000 piece(s)/bag

Product characteristics

Graduation	no
Cap	standard screw cap
Closure type	screw cap
Base shape	flat base
Volume of work	7 ml
Centrifugation max (RCF)	4000 x g

Size

Diameter	13 mm
Length of product	82 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	1000
Piece(s) / pallet	56000
Case volume	0.029 cbm
Weight of product	0.004 kg
Weight of case	4.49 kg
EAN of inner box	4038917082419
EAN case	4038917032513

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

