



Product description

| | |
|---------------------|---|
| Order number | 60.610.001 |
| Product description | Screw cap tube, working volume: 10 ml, (LxØ): 92 x 15.3 mm, material: PP, round base, transparent, screw cap, transparent, cap assembled, sterile, 100 piece(s)/bag |

Product characteristics

| | |
|--------------------------|------------|
| Graduation | no |
| Closure type | screw cap |
| Base shape | round base |
| Volume of work | 10 ml |
| Centrifugation max (RCF) | 3000 x g |

Size

| | |
|-------------------|---------|
| Diameter | 15.3 mm |
| Length of product | 92 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | transparent |
| Tube material | Polypropylene (PP) |

Purity & certification

| | |
|---------------------------|-------------------------|
| Satisfies the requirement | IATA |
| Product category | no medical device IVD |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 22000 |
| Depth of case | 568 mm |
| Width of case | 348 mm |
| Height of case | 155 mm |
| Case volume | 0.0306 cbm |
| Weight of product | 0.0052 kg |
| Weight of case | 3.12 kg |
| EAN of inner box | 4038917114639 |
| EAN case | 4038917014335 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

