



### Product description

Order number	60.610.023
Product description	Screw cap tube, working volume: 10 ml, (LxØ): 92 x 15 mm, material: PP, round base, transparent, screw cap, without cap, 100 piece(s)/bag

### Product characteristics

Graduation	no
Closure type	screw cap
Base shape	round base
Volume of work	10 ml
Centrifugation max (RCF)	4000 x g
Autoclavable (max. °C)	121 °C

### Size

Diameter	15 mm
Height of product	92 mm
Length of product	92 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	50000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	30000
Case volume	0.046 cbm
Weight of product	0.0038 kg
Weight of case	4.26 kg
EAN of inner box	4038917234160
EAN case	4038917234153

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

