



### **Product description**

Order number	60.9922.115
Product description	Screw cap tube, working volume: 25 ml, (LxØ): 54 x 27 mm, material: PS, flat base, transparent, screw cap, yellow, cap assembled, sterile, 50 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	screw cap
Base shape	flat base
Volume of work	25 ml

### **Size**

Diameter	27 mm
Length of product	54 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow

**Purity & certification**

Satisfies the requirement	ADR
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	14000
Depth of case	458 mm
Width of case	373 mm
Height of case	238 mm
Case volume	0.0407 cbm
Weight of product	0.0082 kg
Weight of case	4.8 kg
EAN of inner box	4038917151238
EAN case	4038917027731

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

