



Product description

| | |
|---------------------|---|
| Order number | 62.470 |
| Product description | Tube, working volume: 12 ml, (LxØ): 95 x 16.5 mm, material: PS, flat base, transparent, push cap, without cap, with print, label/print: white, with graduations, 500 piece(s)/bag |

Product characteristics

| | |
|--------------------------|------------|
| Label/ Print | with print |
| Colour of print/label | white |
| Graduation | yes |
| Closure type | push cap |
| Base shape | flat base |
| Volume of work | 12 ml |
| Centrifugation max (RCF) | 3500 x g |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|-------------------|---------|
| Diameter | 16.5 mm |
| Length of product | 95 mm |

Material & colours

| | |
|-------------------|------------------|
| Product material | Polystyrene (PS) |
| Colour of product | transparent |

Purity & certification

| | |
|---------------------------|--------------------------------|
| Satisfies the requirement | sealed at 0.6 bar/10 min |
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 36000 |
| Depth of case | 398 mm |
| Width of case | 398 mm |
| Height of case | 300 mm |
| Case volume | 0.0475 cbm |
| Weight of product | 0.0055 kg |
| Weight of case | 6.21 kg |
| EAN of inner box | 4038917081788 |
| EAN case | 4038917000895 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

