



### Product description

Order number	62.476.322
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PS, round base, transparent, push cap, with print, label/print: white, with graduations, 1,000 piece(s)/bag

### Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	3500 x g

### Size

Diameter	12 mm
Height of product	75 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	80000
Case volume	0.0382 cbm
Weight of product	0.0028 kg
Weight of case	6.14 kg
EAN of inner box	4038917368032
EAN case	4038917368025

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

