




Product Specification
Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, cap enclosed, 1,000 piece(s)/bag

Page 1



Product description

Order number	62.504.037
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 75 x 13 mm, round base, transparent, material: PP, with print, label/print: blue, with graduations, cap enclosed, 1,000 piece(s)/bag, 1,000 piece(s)/carton

Product characteristics

Label/ Print	with print
Colour of print/label	blue
Graduation	yes
Closure type	screw cap
Base shape	round base
Volume of work	5 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round
base, PP, cap enclosed, 1,000 piece(s)/bag**

Size

Diameter	13 mm
Length of product	75 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Tube material	Polypropylene (PP)

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	1000
Piece(s) / pallet	56000
Depth of case	568 mm
Width of case	433 mm
Height of case	120 mm
Case volume	0.0295 cbm
Weight of product	0.0038 kg
Weight of case	4.39 kg
EAN of inner box	4038917412940
EAN case	4038917412933

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

