



### Product description

|                     |   |
|---------------------|---|
| Order number        | 62.526.028  |
| Product description | Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, transparent, push cap, natural, cap assembled, with print, label/print: white, with graduations, sterile, 1 piece(s)/blister |

### Product characteristics

|                          |                 |
|--------------------------|-----------------|
| Label/ Print             | with print      |
| Colour of print/label    | white           |
| Graduation               | yes             |
| Cap                      | Ventilation cap |
| Closure type             | push cap        |
| Base shape               | round base      |
| Volume of work           | 5 ml            |
| Centrifugation max (RCF) | 5000 x g        |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Size**

|                   |       |
|-------------------|-------|
| Diameter          | 12 mm |
| Length of product | 75 mm |

**Material & colours**

|                   |                                  |
|-------------------|----------------------------------|
| Product material  | Polypropylene (PP)               |
| Colour of product | transparent                      |
| Closure material  | Low Density Polyethylene (LD-PE) |
| Colour of cap     | natural                          |

**Purity & certification**

|                           |                          |
|---------------------------|--------------------------|
| Satisfies the requirement | sealed at 0.6 bar/10 min |
| Product category          | no medical device   IVD  |
| Purity standard           | sterile                  |
| Sterilisation             | Electron irradiation     |
| Batched                   | yes                      |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 500           |
| Type of smallest subpackaging | blister       |
| Piece(s) / inner box          | 100           |
| Piece (s) / outer case        | 500           |
| Piece(s) / pallet             | 18000         |
| Depth of case                 | 578 mm        |
| Width of case                 | 368 mm        |
| Height of case                | 195 mm        |
| Case volume                   | 0.0415 cbm    |
| Weight of product             | 0.0029 kg     |
| Weight of case                | 2.82 kg       |
| EAN of inner box              | 4038917109475 |
| EAN case                      | 4038917012492 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

