




**Product Specification**  
**Screw cap tube, 13 ml, (LxØ): 101 x 16.5 mm, PP,  
with print**

**Page 1**



**Product description**

Order number	62.541.062
Product description	Screw cap tube, working volume: 13 ml, (LxØ): 101 x 16.5 mm, material: PP, round base, transparent, screw cap, natural, cap enclosed, with print, label/print: white, with graduations, 500 piece(s)/bag

**Product characteristics**

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	screw cap
Base shape	round base
Volume of work	13 ml
Centrifugation max (RCF)	4000 x g
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 13 ml, (LxØ): 101 x 16.5 mm, PP,  
with print**

**Size**

Diameter	16.5 mm
Length of product	101 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

**Purity & certification**

Satisfies the requirement	IATA
Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	28000
Case volume	0.065 cbm
Weight of product	0.0067 kg
Weight of case	7.2 kg
EAN of inner box	4038917110693
EAN case	4038917534291

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

