



**Product Specification**  
**Screw cap tube, 10 ml, (LxØ): 79 x 16 mm, PP,  
with print**

**Page 1**



**Product description**

Order number	62.551.022
Product description	Screw cap tube, working volume: 10 ml, (LxØ): 79 x 16 mm, material: PP, round base with skirted base, brown, screw cap, brown, cap enclosed, with print, label/print: white, with graduations, 500 piece(s)/bag

**Product characteristics**

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	screw cap
Base shape	round base with skirted base
Volume of work	10 ml
Centrifugation max (RCF)	4000 x g
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Size**

Diameter	16 mm
Length of product	79 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	brown
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown

**Purity & certification**

Satisfies the requirement	IATA
Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	30000
Case volume	0.046 cbm
Weight of product	0.0054 kg
Weight of case	5.9 kg
EAN of inner box	4038917112819
EAN case	4038917030816

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

