



**Product Specification**  
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP,**  
**with print**

**Page 1**



**Product description**

Order number	62.554.100
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PP, conical base, transparent, screw cap, red, cap assembled, with print, label/print: white/blue, with graduations, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 50 piece(s)/rack

**Product characteristics**

Label/ Print	with print
Colour of print/label	white/blue
Graduation	yes
Cap	standard screw cap
Closure type	screw cap
Base shape	conical base
Volume of work	15 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP,  
with print**

**Page 2**

**Size**

Diameter	17 mm
Length of product	120 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	red

**Purity & certification**

Satisfies the requirement	IATA
Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

**Packaging**

Minimum order qty.	25000
Type of smallest subpackaging	rack
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	18000
Case volume	0.0436 cbm
Weight of product	0.0074 kg
Weight of case	4.56 kg
EAN of inner box	4038917376341
EAN case	4038917376334

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

