



**Product Specification**  
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP,  
with print**

**Page 1**



**Product description**

|                     |  |
|---------------------|--|
| Order number        | 62.554.101   |
| Product description | Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PP, conical base, transparent, screw cap, red, cap enclosed, with print, label/print: white/blue, with graduations, 50 piece(s)/bag |

**Product characteristics**

|                          |              |
|--------------------------|--------------|
| Label/ Print             | with print   |
| Colour of print/label    | white/blue   |
| Graduation               | yes          |
| Closure type             | screw cap    |
| Base shape               | conical base |
| Volume of work           | 15 ml        |
| Centrifugation max (RCF) | 20000 x g    |
| Autoclavable (max. °C)   | 121 °C       |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP,  
with print**

**Page 2**

**Size**

|                   |        |
|-------------------|--------|
| Diameter          | 17 mm  |
| Length of product | 120 mm |

**Material & colours**

|                   |                                   |
|-------------------|-----------------------------------|
| Product material  | Polypropylene (PP)                |
| Colour of product | transparent                       |
| Closure material  | High Density Polyethylene (HD-PE) |
| Colour of cap     | red                               |

**Purity & certification**

|                           |                         |
|---------------------------|-------------------------|
| Satisfies the requirement | IATA                    |
| Product category          | no medical device   IVD |
| Batched                   | yes                     |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 500           |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 50            |
| Piece (s) / outer case        | 500           |
| Piece(s) / pallet             | 20000         |
| Case volume                   | 0.0382 cbm    |
| Weight of product             | 0.007 kg      |
| Weight of case                | 3.99 kg       |
| EAN of inner box              | 4038917376310 |
| EAN case                      | 4038917376303 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-21

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**