



Product Specification
Multi-purpose container, 120 ml, (LxØ): 105 x 44 mm, graduated, PP, transparent

Page 1



Product description

| | |
|---------------------|--|
| Order number | 75.9922.420 |
| Product description | Multi-purpose container, max. working volume: 120 ml, (LxØ): 105 x 44 mm, transparent, cap: yellow, graduated, material: PP, screw cap, cap assembled, sterile, 250 piece(s)/bag |

Product characteristics

| | |
|----------------|-----------|
| Graduation | yes |
| Closure type | screw cap |
| Base shape | flat base |
| Volume of work | 120 ml |

Size

| | |
|-------------------|--------|
| Sample volume | 120 ml |
| Diameter | 44 mm |
| Height of product | 105 mm |
| Length of product | 105 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | yellow |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Purity standard | sterile |
| Sterilisation | Gamma irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 250 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 250 |
| Piece (s) / outer case | 250 |
| Piece(s) / pallet | 4000 |
| Depth of case | 548 mm |
| Width of case | 408 mm |
| Height of case | 373 mm |
| Case volume | 0.0834 cbm |
| Weight of product | 0.0188 kg |
| Weight of case | 5.62 kg |
| EAN of inner box | 4038917151863 |
| EAN case | 4038917038461 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

