



Product Specification
Multi-purpose container, 250 ml, (LxØ): 78 x 70 mm, graduated, PP, transparent

Page 1



Product description

| | |
|---------------------|---|
| Order number | 75.9922.532 |
| Product description | Multi-purpose container, max. working volume: 250 ml, (LxØ): 78 x 70 mm, transparent, cap: transparent, graduated, material: PP, screw cap, cap assembled, 240 piece(s)/bag |

Product characteristics

| | |
|----------------|-----------|
| Graduation | yes |
| Closure type | screw cap |
| Base shape | flat base |
| Volume of work | 250 ml |

Size

| | |
|-------------------|--------|
| Sample volume | 250 ml |
| Diameter | 70 mm |
| Height of product | 78 mm |
| Length of product | 78 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Multi-purpose container, 250 ml, (LxØ): 78 x 70 mm, graduated, PP, transparent

Page 2

Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | transparent |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 240 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 240 |
| Piece (s) / outer case | 240 |
| Piece(s) / pallet | 1440 |
| Depth of case | 646 mm |
| Width of case | 450 mm |
| Height of case | 498 mm |
| Case volume | 0.1448 cbm |
| Weight of product | 0.0404 kg |
| Weight of case | 10.93 kg |
| EAN of inner box | 4038917152105 |
| EAN case | 4038917350891 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-12

This document was prepared by EDP support and is valid without signature.



SARSTEDT