



### Product description

Order number	65.3716
Product description	Screw cap, natural, suitable for screw cap micro tubes, for automated filling of screw cap micro tubes and direct adaptation to analyzers, 500 piece(s)/bag

### Product characteristics

Cap	membrane cap
Closure type	screw cap
Suitable for	screw cap micro tubes

### Material & colours

Product material	Polypropylene (PP)
Colour of product	natural
O-ring material	Ethylene propylene diene rubber (EPDM)
Membrane material	rubber

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	648000
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.0004 kg
Weight of case	0.95 kg
EAN of inner box	4038917160193
EAN case	4038917160186

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

