



**Product Specification**  
**Screw cap, violet, suitable for screw cap micro tubes**

**Page 1**



**Product description**

Order number	65.716.008
Product description	Screw cap, violet, suitable for screw cap micro tubes, 500 piece(s)/bag

**Product characteristics**

Cap	flat cap
Closure type	screw cap
Suitable for	screw cap micro tubes

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	violet
O-ring material	Ethylene propylene diene rubber (EPDM)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap, violet, suitable for screw cap micro tubes**

**Page 2**

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	10000
Piece(s) / pallet	660000
Depth of case	438 mm
Width of case	268 mm
Height of case	155 mm
Case volume	0.0182 cbm
Weight of product	0.0003 kg
Weight of case	3.52 kg
EAN of inner box	4038917118613
EAN case	4038917016261

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-07

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**