



**Product Specification**  
**Screw cap, pink, sterile, suitable for screw cap micro tubes**

**Page 1**



**Product description**

Order number	65.716.324
Product description	Screw cap, pink, sterile, suitable for screw cap micro tubes, 500 piece(s)/double bag

**Product characteristics**

Cap	flat cap
Closure type	screw cap
Suitable for	screw cap micro tubes

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	pink
O-ring material	Ethylene propylene diene rubber (EPDM)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Screw cap, pink, sterile, suitable for screw cap  
micro tubes**

**Page 2**

**Purity & certification**

Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Type of smallest subpackaging	double bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	540000
Depth of case	351 mm
Width of case	280 mm
Height of case	138 mm
Case volume	0.0136 cbm
Weight of case	2.05 kg
EAN of inner box	4038917923446
EAN case	4038917923453

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**