



Product Specification
Screw cap, blue, sterile, suitable for screw cap micro tubes

Page 1



Product description

Order number	65.716.723
Product description	Screw cap, blue, sterile, suitable for screw cap micro tubes, 500 piece(s)/double bag

Product characteristics

Cap	flat cap
Closure type	screw cap
Suitable for	screw cap micro tubes

Material & colours

Product material	Polypropylene (PP)
Colour of product	blue
O-ring material	Ethylene propylene diene rubber (EPDM)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**Screw cap, blue, sterile, suitable for screw cap
micro tubes**

Page 2

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	5000
Type of smallest subpackaging	double bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	390000
Depth of case	418 mm
Width of case	278 mm
Height of case	130 mm
Case volume	0.0151 cbm
Weight of product	0.0003 kg
Weight of case	2.01 kg
EAN of inner box	4038917118972
EAN case	4038917053600

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-14

This document was prepared by EDP
support and is valid without signature.



SARSTEDT